

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-NAP1L2-309

Anti NAP1L2

Product type Primary antibodies

Immunogen NAP1L2 (AAH26325.1, 1 a.a. ~ 460 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAKSENKELSESSQEEAGNQIMVEGLGEHLERGEDAAAGLGDDGKCGEEAAAGLGEEGE
NGEDTAAGSGEDGKKGGDTDEDSEADRPKGLIGYVLDTDFVESLPVKVKYRVLALKKLQT
RAANLESKFLREFHDIERKFAEMYQPLLEKRRQIINAIYEPTEEECYKSDSEDCDDEEM
CHEEMYGNEEGMVHEYVDEDDGYEDYYYDYAVEEEEEEEEEEDDIEATGEENKEEEDPKGI
PDFWLTVLKNVDLTPLIKKYDEPILKLLTDIKVKLSDPGEPLSFTLEFHKPNEYFKNE
LLTKTYVLKSKLAYYDHPYRGTAIEYSTGCEIDWNEGKNVTLKTIKKKQKHRIWGTIRT
VTEDFPKDSFFNFFSPHGITSNGRDGNDDFLLGHNLRTYIIPRSVLFFSGDALESQQEGV
VREVNDAIYDKIIYDNWMAAIEEVKACCKNLEALVEDIDR

Gene Alias BPX, MGC26243

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.50 kDa using NAP1L2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of NAP1L2 expression in transfected 293T cell line by NAP1L2 rabbit polyclonal antibody.

Lane 1: NAP1L2 transfected lysate(52.50 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NAP1L2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of NAP1L2 transfected lysate using NAP1L2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NAP1L2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619