

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-NFKBIL2-329

Anti NFKBIL2

Product type	Primary antibodies
Immunogen	NFKBIL2 (AAH08782.1, 1 a.a. ~ 1219 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MRTRLYNLGLTFESLQQTALCNDYFRKSIFLAEQNHLYEDLFRARYNLGTIHWRAQHS QAMRCLEGARECAHTMRKRFMESECCVIAQVLQDLGDFLAALKRALKKAYRLGSQKPVQR AAICQNLQHVLA VVR LQQLEEA EGRDPQGAMVICEQLGDLFSKAGDFPRAAEAYQKQLR FAELDRPGAERAIHVSLATTLGDMKDHGAVRH YEEELRLRSGNVLEEAKTWLNIALS REEAGDAYELLAPCFQKALSCAQQAQRPQLQRQVLQHLHTVQLRLQPQEAPETETRLREL SVAEDEDEEEEAEEAAATAESEALEAGEVELSEGEDDDTGLTPQLEEDEELQGH LGRRKG SKWNRNDMGETLLHRACIEGQLRRVQDLVRQGHPLNPRDYCGWTP LHEACNYGHLEIVR FLLDHGAAVDDPGGQGCEGITPLHDALNCGHFEVAELLER GASVTLRTRKGLSPLETQ QWVKLYRRDL DLETRQKARAMEMLLQAAASGQDPHSSQAFHTPSSLLFDPETSPPLSPCP EPPSNSTR LPEASQVHVRVSPGQAAPAMARPRRSRHGPASSSSSSSEGEDSAGPARPSQKR PRCSATAQRVAAWTPGPASNREAATAST SRAAYQAAIRGVGSAQSRLGPGPPRGH SKALA PQAALIP EEECLAGDWLELD MPLTRSR RPRPRGTGDNRRPSSTSGSDSEESRPRARAKQV RLTCMQSCSAPVNAGPSSLA SEPPGSPSTPRVSEPSGDSSAAGQPLGPAPPPPIRVRVQV QDHLFLIPVPHSSDTHSVAWLAEQAAQRY YQTCGLLPRLTLRKEGALLAPQDLIPDVLQS NDEVLA EVTSWDL PPLTD RYRRACQSLGQGEHQV LQAVELQGLGLSFSACSLALDQAQL TPLLRALKLHTALRELRLAGNRLGDKCVAELVAALGTMPSLALLDLSSNHLGPEGLRQLA MGLPGQATLQSLEELDLSMNPLGDGCGQSLASLLHACPLLSTLRLQACGFGPSFFLSHQT ALGSAFQDAEHLKTL SLSYNALGAPALARTLQSLPAGTLLHLELSSVAAGKGDSDLMEPV FRYLAKEGCALAHLTLSANHLGDKAVRDL CRCLSLCPSLISLDLSANPEISCASLEELLS TLQKR PQGLSFLGLSGCAVQG PLGLGLWDKIAAQLRELQLCSRRLCAEDRDALRQLQPSR PGPGECTLDHGSKLFFRRL


Gene Alias FLJ40087, IKBR

**Application
notes**

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~132.90 kDa using NFKBIL2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NFKBIL2 expression in transfected 293T cell line by NFKBIL2 rabbit polyclonal antibody

Lane 1: NFKBIL2 transfected lysate(132.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619