

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-NPLOC4-176

Anti NPLOC4

Product type	Primary antibodies
Immunogen	NPLOC4 (AAH25930.1, 1 a.a. ~ 608 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAESIIIRVQSPDGVKRIATKRETAATFLKKVAKEFGFQNNGFSVYINRNKTGEITASS NKSLNLLKIKHGDLLFLFPSSLAGPSSEMETSVPPGFKVFGAPNVVEDEIDQYLSKQDGK IYRSRDPQLCRHGPLGKCVHCVPLEPFDELYLNHLEPPVKHMSFHAYIRKLTGGADKGKF VALENISCKIKSGCEGHPWPNGICTKCQPSAITLNRQKYRHVDNIMFENHTVADRFLDF WRKTGNQHFGYLYGRYTEHKDIPLGIRAEVAAIYEPPIGTQNSLELLEDPKAEVVDEIA AKLGLRKVGWIFTDLVSEDTRKGTVRYSRNKDTYFLSSEECITAGDFQNKHPNMCRLSPD GHFGSKFVTAVATGGPDNQVHFEGYQVSNQCMALVRDECLLPCKDAPELGYAKESSEQY VPDVFYKDVDFKFGNEITQLARPLPVEYLIIDITTFPKDPVYTFSSISQNPFIENRDVLG ETQDFHSLATYLSQNTSSVFLDTISDFHLLFLVTNEVMPLQDSISLLLEAVRTRNEELA QTWKRSEQWATIEQLCSTVGGQLPGLHEYGAVGGSTHTATAAMWACQHCTFMNQPGTGHC EMCSLPRT
Gene Alias	FLJ20657, FLJ23742, KIAA1499, NPL4

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~68.10 kDa using NPLOC4 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NPLOC4 expression in transfected 293T cell line by NPLOC4 rabbit polyclonal antibody.


Lane 1: NPLOC4 transfected lysate(68.10 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NPLOC4 transfected lysate.

Immunoprecipitation



Immunoprecipitation of NPLOC4 transfected lysate using NPLOC4 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NPLOC4 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619