

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-NSFL1C-189

Anti NSFL1C

Product type	Primary antibodies
Immunogen	NSFL1C (NP_057227.2, 1 a.a. ~ 370 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAERQEALREFVAVTGAEEDRARFFLESAGWDLQIALASFYEDGGDEDIVTISQATPSS VSRGTAPSDNRVTSFRDLIHDQDEDEEEEEEGQRFYAGGSERSGQQIVGPPRKKSPNELVD DLFKGAKEHGAVAVERTKSPGETSKPRPFAGGGYRLGAAPEEESAYVAGEKRQHSSQDV HVVLLKLWKSGFSLDNGELRSYQDPSNAQFLESIRRGEPVPAELRRLAHGGQVNLDMEDHRD EDFVKPKGAFKAFTGEGQKLGSTAPQVLSTSSPAQQAENEAKASSSILIDSEPTTNIQI RLADGGRLVQKFNHSHRISDIRLFIVDARPAMAATSFILMTTFPNKELADESQTLEKANL LNAVIVQRLT
Gene Alias	MGC3347, UBX1, UBXD10, UBXN2C, dJ776F14.1, p47
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.60 kDa using NSFL1C transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NSFL1C expression in transfected 293T cell line by NSFL1C rabbit polyclonal antibody.

Lane 1: NSFL1C transfected lysate(40.60 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619