



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-OIP5-392

Anti OIP5

Product type	Primary antibodies
Immunogen	OIP5 (NP_009211.1, 1 a.a. ~ 229 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MAAQPLRHRSRCATPPRGDFCGGTERAIDQASFTTSMEWDTQVVKGSSPLGPAGLGAEPP
AAGPQLPSWLQPERCAVFQCAQCHAVLADSVHLAWDLRSRLGAVVFSRVTNVNVLEAPFL
VGIEGSLKGSTYNLLFCGSCGIPVGFHLYSTHAALAALRGHFCLSSDKMVCYLLKTKAIV
NASEMDIQNVPLSEKIAELKEKIVLTHNRLKSLMKILSEVTPDQSKPEN


5730547N13Rik, LINT-25, MIS18beta

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~24.70 kDa using OIP5 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of OIP5 expression in transfected 293T cell line by OIP5 rabbit polyclonal antibody.

Lane 1: OIP5 transfected lysate(24.70 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of OIP5 transfected lysate.

Immunoprecipitation



Immunoprecipitation of OIP5 transfected lysate using OIP5 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with OIP5 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619