



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-ORC6L-547

Anti ORC6L

Product type Primary antibodies

Immunogen ORC6L (AAH63565.1, 1 a.a. ~ 252 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MGSELIGRLAPRLGLAEPDMLRKAE EYLRLSRVKCVGLSARTTETSSAVMCLDLAASWMK
CPLDRAYLIKLSGLNKETYQSCLKSFECLLGLNSNIGIRDLAVQFSCIEAVNMASKILKS
YESSLPQTQQVDLDLSRPLFTSAALLSACKILKLV DKNKMMVATSGVKKAI FDR LCKQLE
KIGQQVDREPGDVATPPRKRKKIVVEAPAKEME KVEEMPHKPQKDEDLTQDYEEWKRKIL
ENAAS AQKATAE


Gene Alias ORC6

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~27.83 kDa using ORC6L transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of ORC6L expression in transfected 293T cell line by ORC6L rabbit polyclonal antibody.


Lane 1: ORC6L transfected lysate(27.83 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ORC6L transfected lysate

Immunoprecipitation



Immunoprecipitation of ORC6L transfected lysate using ORC6L rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ORC6L purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619