

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-PAFAH2-101

Anti PAFAH2

Product type	Primary antibodies
Immunogen	PAFAH2 (NP_000428.2, 1 a.a. ~ 392 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGVNQSVGFPPVTGPHLVGCGDVMEGQNLQGSFFRLFYPCQKAEETMEQPLWIPRYEYCT GLAEYLQFNKRCGGLLFNLAVGSCRLPVSWNGPFFKTKDSGYPLIIFSHGLGAFRTLVSF CMELASRGFVVAVPEHRDRSAATTYFCKQAPEENQPTNESLQEEWIPFRRVEEGEKEFHV RNPQVHQRVSECLRVLKILQEVTAGQTVFNILPGGLDMLTKGNIDMSRVAVMGHSFGGA TAILALAKETQFRCAVALDAWMFPLERDFYPKARGPVFFINTEKFQTMESVNLMMKKICAQ HEQSRIITVLGSHVRSQTDFAFVTGNLIGKFFSTETRGS�DPYEGQEVMMVRAMLAFLLKH LDLKEDYNQWNNLIEGIGPSLTPGAPHHSSL
Gene Alias	FLJ26025, HSD-PLA2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.00 kDa using PAFAH2 transfected 293T lysate.

..
Western Blotting (Transfected lysate)



Western Blot analysis of PAFAH2 expression in transfected 293T cell line by PAFAH2 rabbit polyclonal antibody

Lane 1: PAFAH2 transfected lysate(44.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619