

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PAICS-391

Anti PAICS

Product Primary antibodies

type

Immunogen PAICS (NP_006443.1, 1 a.a. ~ 425 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MATAEVLNIGKKLYEGKTKEVYELLDSPGKVLLQSKDQITAGNAARKNHLEGKAAISNKI
TSCIFQLLQEAGIKTAFTKCGETAFIAPQCEMPIEWVCRRIATGSFLKRNPVGVEGYK
FYPPKVELFFKDDANNDPQWSEEQLIAAKFCFAGLLIGQTEVDIMSHATQAIFEILEKSW
LPQNCTLVDMKIEFGVDVTTKEIVLADVIDNDSWRLWPSGDRSQQKDKQSYRDLKEVTPE
GLQMVKKNFVVAERVELLLKSESQCRVVLMGSTDLGHCEKIKKACGNFGIPCELVRT
SAHKGPDETLRIKAEYEGDIPTVFVAVAGRSNGLGPVMSGNTAYPVISCPPLTPDWGVQ
DVWSSLRLPSGLGCSTVLSPEGSAQFAAQIFGLSNHLVWSKLRASILNTWISLKQADKKI
RECNL

ADE2, ADE2H1, AIRC, DKFZp781N1372, MGC1343, MGC5024, PAIS

Application WB

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.10 kDa using PAICS transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of PAICS expression in transfected 293T cell line by PAICS rabbit polyclonal antibody.

Lane 1: PAICS transfected lysate(47.10 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619