



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PCID2-397

Anti PCID2

Product type	Primary antibodies
Immunogen	PCID2 (AAH31246.1, 1 a.a. ~ 397 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MRLTDVVQQLVYEAIIDSRDGASCAELVSFKHPHVANPRLQMASPEEKCCQQVLEPPYDEMFAAHLRCTYAVGDHDFIEAYKCQTVIVQSFLRAFQAHKEENWALPVMYAVALDLRVFANNA
DQQLVKKGKSKVGDMLLEKAAELLMSCFRVCASDTRAGIEDSKKWGMLFLVNQLFKIYFKI
NKLHLCKPLIRAIIDSSNLKDDYSTAQRVTYKYYVGRKAMFDSDFKQAEYLSFAFEHCHR
SSQKNRMILYLLPVKMLLGHMPTVELLKKYHLMQFAEVTRAVSEGNLLLLHEALAKHE
AFFIRCGIFLILEKLIITYRNLFKKVYLLKTHQLSLDAFLVALKFMQVEDVDIDEVQC
ILANLIYMGHVKGYSIQHQKLVVSKQNPFPPLSTVC


DKFZp686C20226, F10, FLJ11305, FLJ99362, MGC16774

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.80 kDa using PCID2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PCID2 expression in transfected 293T cell line by PCID2 rabbit polyclonal antibody.

Lane 1: PCID2 transfected lysate(45.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619