

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-PDE4A-234

Anti PDE4A

Product type	Primary antibodies
Immunogen	PDE4A (NP_006193.1, 1 a.a. ~ 647 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPLVDFFCETCSKPWLVGWWDQFKRMLNRELTHLSEMSRSGNQVSEYISTTFLDKQNEVE IPSPMKEREKQQAPRPRPSQPPPPVPHLQPMSTGLKMLMHSNSLNNSNIPRFGVKT DQEELLAQELNENLKWGLNIFCVSDYAGGRSLTCIMYMIFQERDLLKKFRIPVDTMVTYM LTLEDHYHADVAYHNSLHAADVQLQSTHVLLATPALDAVFTDLEILAALFAAAIHVDVHPG VSNQFLINTNSELALMYNDESVLENHHLAVGFKLLQEDNCDIFQNLKRQRQSLRKMVID MVLATDMSKHMILLADLKTMTVETKKVTSSGVLLLDNYSDRIQVLRNMVHCADLSNPTKPL ELYRQWTDRIAMAEFFQQGDRERERGMEISPMCDKHTASVEKSQVGFIDYIVHPLWETWAD LVHPDAQEILDTELDNRDWYYSAIRQSPSPPEEESRGPGHPLPDKFQFELTLEEEEEEE EISMAQIPCTAQEALTAQGLSGVEEALDATIAWEASPAQESLEVMAQEASLEAELEAVYL TQQAQSTGSAPVAPDEFSSREEFVAVSHSSPSALALQSPLLPAWRTLVSSEHAPGLPGL PSTAAEVEAQREHQAAKRACSACAGTFGEDTSALPAPGGGGSGGDPT
Gene Alias	DPDE2, PDE4, PDE46
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~72.20 kDa using PDE4A transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PDE4A expression in transfected 293T cell line by PDE4A rabbit polyclonal antibody.

Lane 1: PDE4A transfected lysate(72.20 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619