



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-PLA1-187

Anti PLA1A

Product type	Primary antibodies
Immunogen	PLA1A (NP_056984.1, 1 a.a. ~ 456 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPPGPWESCFWVGGLILWLSVGSSGDAPPTPQPKCADFQSANLFEGTDLKVQFLLFVPSN PSCGQLVEGSSDLQNSGFNATLGTKLIHGFRVLGTKPSWIDTFIRTLRATNANVIAVD WIYGSTGVYFSAVKNVIKLSLEISLFLNKLLVLGVSESSIIHGVSLGAHVGGMVGQLFG GQLGQITGLDPAGPEYTRASVEERLDAGDALFVEAIHTDNDNLGIRIPVGHVDYFVNGGQ DQPGCPTFFYAGYSYLICDHMRVHLYISALENSCPLMAFPCASYKAFLAGRCLDCFNPF LLSCPRIGLVEQGGVKIEPLPKEVKVYLLTTSSAPYCMHHSLEFHLKELRNKDTNIEVT FLSSNITSSSKITIPKQQRYGKGIIAHATPQCQINQVKFKFQSSNRVWKKDRTTIIGKFC TALLPVNDREKMOVCLPEPVNLQASVTVSCDLKIACV
Gene Alias	PS-PLA1, PSPLA1
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.70 kDa using PLA1A transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PLA1A expression in transfected 293T cell line by PLA1A rabbit polyclonal antibody.


Lane 1: PLA1A transfected lysate(49.70 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PLA1A transfected lysate.

Immunoprecipitation



Immunoprecipitation of PLA1A transfected lysate using PLA1A rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PLA1A mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619