



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-POLR1D-036

Anti POLR1D


Product type	Primary antibodies
Immunogen	POLR1D (NP_057056.1, 1 a.a. ~ 133 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MEEDQELERKISGLKTSMAEGERKTALEMVQAAGTDRHCVTFVLHEEDHTLGNSLRYMIM
KNPEVEFCGYTTTHPSESINLRIQTRGTLPAVEPFQRLNELMNVQCQHVLDKFEASIKD
YKDQKASRNESTF

Gene Alias FLJ20616, MGC9850, POLR1C, RPA16, RPA9, RPAC2, RPO1-3

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~15.20 kDa using POLR1D transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of POLR1D expression in transfected 293T cell line by POLR1D rabbit polyclonal antibody.


Lane 1: POLR1D transfected lysate(15.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of POLR1D transfected lysate.

Immunoprecipitation



Immunoprecipitation of POLR1D transfected lysate using POLR1D rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with POLR1D rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619