

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-POLR3C-297

Anti POLR3C

Product type	Primary antibodies
Immunogen	POLR3C (NP_006459.3, 1 a.a. ~ 534 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MTQAEIKLCSLLLQEHFGEIVEKIGVHLIRTGSQPLRVIAHDTGTSLDQVKKALCVLVQH NLVSYQVHKRGVVEYEAQCSRVLRLRYPRYIYTTKTLYSDTGELIVEELLNGKLTMSA VVKKVADRLTETMEDGKTMDYAEVSNTFVRLADTHFVQRCPSVPTTENS DPGPPPPAPTL VINEKDMYLVPKLSLIGKGKRRRSSDEDAAGEPKAKRPKYTTDNKEPIPDDGIYWQANLD RFHQHFRDQAIVSAVANRMDQTSSEIVRTMLRMSEITTSSAPFTQPLSSNEIFRSLPVG YNISKQVLDQYLTLADDPLEFVGKSGDSGGGMYVINLHKALASLATATLESVVQERFGS RCARIFRLVLQKKHIEQKQVEDFAMIPAKEAKDMLYKMLSENFMSLQEIPKTPDHAPSRT FYLYTVNILSAARMLLHRCYKSANLIERRQFETKENKRLLEKSQRVEAIIASMQATGAE EAQLQEIEEMITAPERQQLETLKRVNKLDAEIQVDETIFLLESYIECTMKRQ
Gene Alias	RPC3, RPC62
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~60.60 kDa using POLR3C transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of POLR3C expression in transfected 293T cell line by POLR3C rabbit polyclonal antibody.

Lane 1: POLR3C transfected lysate(60.60 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619