

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PPT2-035


Anti PPT2

Product type	Primary antibodies
Immunogen	PPT2 (NP_005146.3, 1 a.a. ~ 302 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MLGLWGQRLPAAWVLLLLPFLPLLLAAPHRASYPVIVVHGLFDSSYSFRHLLLEYIN ETHPGTVVTVDLFDGRESLRPLWEQVQGFREAVVPIMAKAPQGVHLICYSQGGLVCRAL LSVMDDHNVDSFISLSSPQMGGYGDYDLKWLFPSTMRNLYRICYSPWQGQEFSCNYWH DPHHDDLVLNASSFLALINGERDHPNATVWRKNFLRVGHLVLIGGPDDGVITPWQSSFFG FYDANETVLEMEEQLVYLRDSFGLKTLLARGAIVRCMPMAGISHTAWHSNRTLYETCIEPW LS
---------------------------	---

Gene Alias	C6orf8, DKFZp564P1516, G14
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.30 kDa using PPT2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PPT2 expression in transfected 293T cell line by PPT2 rabbit polyclonal antibody.

Lane 1: PPT2 transfected lysate(34.30 kDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619