

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PRMT7-093

Anti PRMT7


Product type	Primary antibodies
Immunogen	PRMT7 (NP_061896.1, 1 a.a. ~ 692 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MKIFCSRANPTTGSVEWLEEDHYHDIHQEIARSSYADMLHDKDRNVKYYQGIRAAVSRVK DRGQKALVLDIGTGTGLLSMMAVTAGADFCYAIEVFKPMADAADVIVEKNGFSDKIKVIN KHSTEVTVGPEGDMPCRANILVTELFDTLIGEGALPSYEHHRHLVEENCEAVPHRATV YAQLVESGRMWSWNKLFPIHVQTSLGQVIVPPVDVESC PGAPSVCDIQLNQVSPADFTV LSDVLPMSIDFSKQVSSSAACHSRRFEPLTSGRAQVVLSSWWDIEMDPEGKIKCTMAPFW AHS DP EEMQWRDHWMQCVYFLPQEEPVVQGSALYLVAHHDDYCVWYSLQRTSPEKNERVR QMRPVCDCQAHLWNRPRFGEINDQDRTDRYVQALRTVLKPDVCLCVSDGSLLSVLAHH LGVEQVFTVESSAASHKLLRKIFKANHLEDKINIIIEKRPELLTNE DLQGRKVSLLLGE PF FTTSLLPWHNLYFWYVRTAVDQHLGPGAMVMPQAASLHVVVEFRDLWRIRSPCGDCEGF DVHIMDDMIKRALDFRESREAPHPLWEYPCRSLSPEWQILTDFDQQPVPLQLCAEGTV ELRRPGQSHA AVLWMEYH LTPECTLSTGLLEPADPEGGCCWNPHCKQAVYFFSPAPDPRA LLGGPRTVSYAVEFHPDTGDIIMEFRHADTPD
---------------------------	--

Gene Alias	FLJ10640, KIAA1933
-------------------	--------------------

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~78.50 kDa using PRMT7 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)




Western Blot analysis of PRMT7 expression in transfected 293T cell line by PRMT7 rabbit polyclonal antibody.

Lane 1: PRMT7 transfected lysate(78.50 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PRMT7 expression in human kidney (Predicted Mw: 78.4590 kDa).

Western Blotting (Tissue lysate)




PRMT7 rabbit polyclonal antibody. Western Blot analysis of PRMT7 expression in human kidney.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PRMT7 transfected lysate.

Immunoprecipitation



Immunoprecipitation of PRMT7 transfected lysate using PRMT7 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PRMT7 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619