

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PSMB4-384

Anti PSMB4

Product type	Primary antibodies
Immunogen	PSMB4 (NP_002787.2, 1 a.a. ~ 264 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEAFLGSRSGLWAGGPAPGQFYRIPSTPDSFMDPASALYRGPIRTQNPMVTGTSVLGVK FEGGVVIAADMLGSYGSLARFRNISRIMRVNNSTMLGASGDYADFQYLKQVLGQMVIDEE LLGDGHSYSPRAIHSWLTRAMYSRRSKMNPLWNTMVIGGYADGESFLGYVDM LGVAYEAP SLATGYGAYLAQPLLREVLEKQPVLSQTEARDLVERCMRVLYYRDARSYNRFQIATVTEK GVEIEGPLSTETNWDIAHMISGFE


HN3, HsN3, PROS26

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~29.20 kDa using PSMB4 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of PSMB4 expression in transfected 293T cell line by PSMB4 rabbit polyclonal antibody.

Lane 1: PSMB4 transfected lysate(29.20 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PSMB4 expression in Jurkat (Predicted Mw: 29.2040 kDa).


Western Blotting (Cell lysate)



PSMB4 rabbit polyclonal antibody. Western Blot analysis of PSMB4 expression in Jurkat.

Antibody dilution: 1:1000. Western Blot analysis of PSMB4 expression in human placenta (Predicted Mw: 29.2040 kDa).

Western Blotting (Tissue lysate)



PSMB4 rabbit polyclonal antibody. Western Blot analysis of PSMB4 expression in human placenta.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619