

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-PSRC1-321

Anti PSRC1

Product type	Primary antibodies
Immunogen	PSRC1 (NP_116025.1, 1 a.a. ~ 333 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEDLEEDVRFIVDETLDGGLSPSDSREEEDITVLVTPEKPLRRGLSHRSDPNAVAPAPQ GVRLSLGPLSPEKLEEILDEANRLAAQLEQCALQDRESAGEGLGPRRVKPSPRRETFLVK DSPVRDLLPTVNSLTRSTPSPSSLTPRLRSNDRKGSVRALRATSGKRPSNMKRESPTCNL FPASKSPASSPLTRSTPPVRGRAGPSGRAAASPPTPIRSVLAPQPSTSNSQRLPRPQGAA AKSSSQLPIPSAIPRPASRMPLTSRSVPPGRGALPPDSLSTRKGLPRPSTAGHRVRESGH KVPVSQRLNLPVMGATRSNLQPPRKVAVPGPTR
Gene Alias	DDA3, FP3214, MGC1780, RP11-297O4.2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.60 kDa using PSRC1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PSRC1 expression in transfected 293T cell line by PSRC1 rabbit polyclonal antibody.

Lane 1: PSRC1 transfected lysate(35.60 kDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619