

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PUS1-304

Anti PUS1

Product type Primary antibodies

Immunogen PUS1 (AAH02901.1, 1 a.a. ~ 400 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAGNAEPPPPAGAACPQDRRSCSGRAGGDRVWEDGEHPAKKLKSGGDEERREKPPKRKIVL
LMAYSGKGYHGMQRNVGSSQFKTIEDDLVSALVRSGCIPENHGDMRKMSFQRCARTDKG
VSAAGQVVSLKVWLIDDILEKINSHLPISHIRILGLKRVTGGFNSKNRCDARTYCYLLPTF
AFAHKDRDVQDETYRLSAETLQQVNRLACYGKTHNFHNFTSQKGPQDPSACRYILEMYC
EFPFVREGLEFAVIRVKGQSFMMHQIRKMVGLVVAIVKGYAPESVLSRWSGTEKVDVPKA
PGLGLVLERVHFEKYNQRFNGDGLHEPLDWAQEEGKVAAFKEEHYPTIIGTERDERSMA
QWLSTLPIHNFSATALTAGGTGAKVPSPLEGSEGDGDTD*

Gene Alias MGC11268, MLASA

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.00 kDa using PUS1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PUS1 expression in transfected 293T cell line by PUS1 rabbit polyclonal antibody.


Lane 1: PUS1 transfected lysate(44.00 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PUS1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of PUS1 transfected lysate using PUS1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PUS1 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619