

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-PWP2-402


Anti PWP2

Product type	Primary antibodies
Immunogen	PWP2 (NP_005040.2, 1 a.a. ~ 919 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKFAYRFSNLLGTVYRRGNLNFCDGNSVISPVGNRVTVFDLKNNKSDTLPLATRYNVKC VGLSPDGR LAIVDEGGDALLVSLVCRSVLHHFHFKGSVHSVSVFSPDGRKFVVTKG NIAQ MYHAPGKKREFNAFVLDKTYFGPYDETT CIDWTDDSRCFVVGSKDMSTWVFGAERWDNLI YYALGGHKDAIVACFFESNSLDLYSLSQDGVLCMWQC DTPPEGLRLKPPAGWKADLLQRE EEEEEEEDQEGDRETTIRGKATPAEEEEKTGKV KYSRLAKYFFNKEGDFNNLTAAAFHKKS HLLVTGFASGIFHLHELPEFNLIHLSISDQSIASVAINSSGDWIAFGCSGLGQLLVWEW QSESYVLKQQGHFNSMVALAYSPDGQYIVTGGDDGKV KVVWNTLSGFCFVTFTEHSSGV TG VTFTATGYVVVTSSMDGTVRAFDLHRYRNFRTFTSPRPTQFSCVAVDASGEIVSAGA QDS FEIFVWSMQTGRLLDVLSGHEGPISGLCFNPMKSVLASASWDKTVRLWDMFDSWRTKETL ALTSDALAVTFRPDGAELAVATLNSQITFWDPENAVQTGSIEGRHDLKTGRKELDKITAK HAAKGKAFTALCYSADGHSILAGGMSKFVCIYHVREQILMKRFEISCNLSLDAMEEFLNR RKMTEFGNLALIDQDAGQEDGVAIPLPGVRKGDMS SRHFKEIRVTS LRFSPTGRCWAAT TTEGLLIYSLDTRVLDFDPFELDTSVTPGRVREALRQQDFTRAILMALRLNESKLVQEAL E AVPRGEIEVVTSSLPELYVEKVLFLASSFEVSRHLEFYLLWTHKLLMLHGQKLKSRAGT LLPVIQFLQKSIQRHLDDLSKLCSWNHYNMQYALAVSKQRGTKRSLDPLGSEEEAEASED DSLHLLGGGGRDSEEEMLA

EHOC-17, PWP2H

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~102.50 kDa using PWP2 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of PWP2 expression in transfected 293T cell line by PWP2 rabbit polyclonal antibody.


Lane 1: PWP2 transfected lysate(102.50 kDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PWP2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of PWP2 transfected lysate using PWP2 rabbit polyclonal antibody and Protein A Magnetic Bead, and Immunoblotted with PWP2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619