

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RAD23B-505

Anti RAD23B

Product type Primary antibodies

Immunogen RAD23B (AAH20973.1, 1 a.a. ~ 409 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MQVTLKTLQQQTFKIDIDPEETVKALKEKIESEKGDAPFPVAGQKLIYAGKILNDDTALK
EYKIDEKNFVVVMVTKPKAVSTPAPATTQQSAPASTTAVTSSTTTTVAQAPTPVPALAPT
STPASITPASATASSEPPAPASAAKQEKPAEKPAETPVATSPTATDSTSGDSSRSNLFEDA
TSALVTGQSYENMVTEIMSMGYEREQVIAALRASFNPDRAVEYLLMGIPGDRESQAVVD
PPQAASTGVPQSSAVAAAAATTTATTTTTSSGGHPLEFLRNQPQFQQMRQIIQQNPSSLP
ALLQQIGRENPLLQQISQHQEHFIQMLNEPVQEAGGQGGGGGGSGGIAEAGSGHMNYI
QVTPQEKEAIERLKALGFPEGLVIQAYFACEKNENLAANFLLQQNFDED


Gene Alias HHR23B, HR23B, P58

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.99 kDa using RAD23B transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RAD23B expression in transfected 293T cell line by RAD23B rabbit polyclonal antibody.

Lane 1: RAD23B transfected lysate(44.99 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619