

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RNF111-338

Anti RNF111

Product type Primary antibodies

Immunogen RNF111 (AAH60862.1, 1 a.a. ~ 985 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MSQWTPEYKELYTLKVDKSEIPSDAPKTQESLKGILLHPEPIGAAKSFPAGVEMINSKV
GNEFSHLCDDSDQKQEKEMNGNQEQEKSLLVVRKKRKSQQAGPSYVQNCVKENQGILGLRQ
HLGTPSDEDNDSSFSDCSSPSSSLHFGSDTVTSDEDKEVSVRHSQTILNAKSRSHSAR
SHKWPRTESETSVSGLLMKRPCLHGSSLRRLPCRKFVKNNSSQRTQKQKERILMQRKKRE
VLARRKYALLPSSSSSSSENDLSSESSSSSSSTEGEEDLFVSAENHQNNPAVPSGSDIEDV
VVIEASSTPQVTANEEINVTSTDSEVEIVTVGESYRSRSTLGHSRSHWSQGSSSHASRPQ
EPRNRSRISTVIQPLRQNAAEVVDLTVDDEPTVPTTSARMESQATSASINNSNPSTSE
QASDTASAVTSSQPSTVSETSATLTSTSTTGTSIGDDSRRTTSSAVTETGPPAMPRLPSC
CPQHSPCGGSSQNHHALGHPHTSCFQQHGHFQHGHGHGHHTPHPAVPVSPSFSDPACPVE
RPPQVQAPCGANSSSGTSYHEQALPVDLSNSGIRSHGSGSFGASAFDPCCPVSSSRAAI
FGHQAAAAAPSQPLSSIDGYGSSMVAQPQPQPPQPSLSSCRHYMPPPYASLTRPLHHQA
SACPHSHGNPPPTQPPPQVDYVIPHPVHAFHSQISSHATSHPVAPPPPTHLASTAAPIP
QHLPPTHQPISSHIPATAPPAQRLHPHEVMQRMVQRRRMMQHPTRAHERPPPHPHRMHP
NYGHGHHIHPVQTMSSHPRQAPERSAWELGIEAGVTAATYTPGALHPLAHYHAPPRLHH
LQLGALPLMVPDMAGYPHIYISSGLDGTSGFRGPFGRNFEELIHLEERLGNVNRGASQGT
IERCTYPHKYKKRKLHCKQDGEEGTEEDTEEKCTICLSILEEGEDVRRLLPCMHLFHQVCV
DQWLITNKKCPICRVDEIAQLPSES

Gene Alias ARK, DKFZp313E0731, DKFZp686H1966, DKFZp761D081, FLJ38008

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~107.80 kDa using RNF111 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of RNF111 expression in transfected 293T cell line by RNF111 rabbit polyclonal antibody.

Lane 1: RNF111 transfected lysate(107.80 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of RNF111 expression in NIH/3T3 (Predicted Mw: 107.7820 kDa)

Western Blotting (Cell lysate)



RNF111 rabbit polyclonal antibody. Western Blot analysis of RNF111 expression in NIH/3T3.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619