

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RPA4-551

Anti RPA4


Product type	Primary antibodies
Immunogen	RPA4 (NP_037479.1, 1 a.a. ~ 261 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSKSGFGSYGSISAADGASGGSDQLCERDATPAIKTQRPKVRIQDVVPCNVNQLLSSTVF DPVFKVRGIIVSQVSIVGVIRGAEKASNHICYKIDDMTAKPIEARQWFGREKVKQVTPLS VGYYVKVFGILKCPTGTKSLEVLKIHVLEDMNEFTVHILETVNAHMMMLDKARRDTTVESV PVSPSEVNDAGDNDESHRNFIQDEVLRLLIHECPHQEGKSIHELRAQLCDLSVKAIKEAID YLTVEGHIYPTVDREHFKSAD

Gene Alias HSU24186, MGC120333, MGC120334

Application **WB**

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~28.90 kDa using RPA4 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of RPA4 expression in transfected 293T cell line by RPA4 rabbit polyclonal antibody.

Lane 1: RPA4 transfected lysate(28.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of RPA4 expression in human kidney (Predicted Mw: 28.8680 kDa).

Western Blotting (Tissue lysate)




RPA4 rabbit polyclonal antibody. Western Blot analysis of RPA4 expression in human kidney.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RPA4 transfected lysate

Immunoprecipitation



Immunoprecipitation of RPA4 transfected lysate using RPA4 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RPA4 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619