

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-SFTPB-593

Anti SFTPB

Product type	Primary antibodies
Immunogen	SFTPB (NP_000533.2, 1 a.a. ~ 381 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAESHLQWLLLLPTLCGPGTAAWTTSSLACAQGPEFWCQSLEQALQCRAIGHCLQEWW GHVGADDLCQECEDIVHILNKMAKEAIFQDTMRKFLEQECNVLPKLLMPQCNQVLDDYF PLVIDYFQNQTDSENGICMHLGLCKSRQPEPEQEPGMSDPLPKPLRDPLPDLLDKLVLPV LPGALQARPGPHTQDLSEQQFPIPLPYCWLCRALIKRIQAMIPKGALAVAVAQVCRVVPL VAGGICQCLAERYSVILLDTLLGRMLPQLVCRLVLRCSMDDSAGPRSPTEGEWLPRDSECH LCMSVTTQAGNSSEQAIPQAMLQACVGSWLDREKCKQFVEQHTPQLLLTVPRGWDAAHTTC QALGVCMTMSSPLQCIHSPDL
Gene Alias	PSP-B, SFTB3, SFTP3, SMDP1, SP-B
Application notes	WB Antibody dilution: 1:1000. Western Blot analysis of SFTPB expression in human pancreas (Predicted Mw: 41.7280 kDa).

Western Blotting (Tissue lysate)




SFTPB rabbit polyclonal antibody. Western Blot analysis of SFTPB expression in human pancreas.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of SFTPb transfected lysate

Immunoprecipitation



Immunoprecipitation of SFTPb transfected lysate using SFTPb rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with SFTPb mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619