

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-SIRT5-087


Anti SIRT5

Product type	Primary antibodies
Immunogen	SIRT5 (NP_036373.1, 1 a.a. ~ 310 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MRPLQIVPSRLISQLYCGLKPPASTRNQICLKMARPSSSMADFRKFFAKAKHIVIISGAG VSAESGVPTFRGAGGYWRKWQAQDLATPLAFAHNPSRVWEFYHYRREVMGSKEPNAGHRA IAECETRLGKQGRRVVITQNIDELHRKAGTKNLLEIHGSLFKTRCTSCGVVAENYKSPI CPALSGKGAPEPGTQDASIPVEKLPRCEEAGCGLLRPHVWVWFGENLDPAILLEEVDRELA HCDLCLVVGTSVVYPAMFAPQVAARGVPVAEFNTETTPATNRFRHFHQGPCGTTLPEA LACHENETVS

Gene Alias FLJ36950, SIR2L5

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~33.90 kDa using SIRT5 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of SIRT5 expression in transfected 293T cell line by SIRT5 rabbit polyclonal antibody.

Lane 1: SIRT5 transfected lysate(33.90 kDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of SIRT5 expression in HeLa (Predicted Mw: 33.8810 kDa).


Western Blotting (Cell lysate)



SIRT5 rabbit polyclonal antibody. Western Blot analysis of SIRT5 expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of SIRT5 expression in human kidney (Predicted Mw: 33.8810 kDa).

Western Blotting (Tissue lysate)



SIRT5 rabbit polyclonal antibody. Western Blot analysis of SIRT5 expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619