

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-SNX9-356

Anti SNX9


Product type	Primary antibodies
Immunogen	SNX9 (NP_057308.1, 1 a.a. ~ 595 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MATKARVMYDFAAEPGNNELTVNEGEIITITNPDVGGGWLEGRNIKGERGLVPTDYVEIL PSDGKQDFSCGNSVADQAFDLSASTAQASSAASNNHQVSGNDPWSAWSASKSGNWE SSEGWGAQPEGAGAQRNTNTPNNWDATFGHPQAYQGPATGDDDDWDEDWDGPKSSSYFKD SESADAGGAQRGNSRASSSSMKIPLNKFPGFAPGTEQYLLAKQLAKPKEKIPPIVGDYG PMWVYPTSTFDCVVADPRKGSKMYGLKSYIEYQLTPTNTNRSVNHRYKHFDWLYERLLVK FGSAIPIPSLPDKQVTGRFEEEFIKMRMERLQAWMTRMCRHPVISESEVFQQFLNFRDEK EWKTGKRKAERDELAVMIFSTMEPEAPDLDLVEIEQKCEAVGKFTKAMDDGVKELLTVG QEHWKRCTGPLPKEYQKIGKALQSLATVFSSSGYQGETDLNDAITEAGKTYEEIASLVAE QPKKDLHFLMECNHEYKGLGCFPDIIIGTHKGAIEKVKESDKLVATSKITLQDKQNMVVR VSIMSYALQAEMNHFHSNRIYDYNVIRLYLEQQVQFYETIAEKLRLQALSFRFPVM
	MST155, MSTP155, SDP1, SH3PX1, SH3PXD3A, WISP

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~66.60 kDa using SNX9 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SNX9 expression in transfected 293T cell line by SNX9 rabbit polyclonal antibody.


Lane 1: SNX9 transfected lysate(66.60 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of SNX9 transfected lysate.

Immunoprecipitation



Immunoprecipitation of SNX9 transfected lysate using SNX9 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with SNX9 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619