

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TBCE-385

Anti TBCE

Product Primary antibodies

type

Immunogen TBCE (NP_003184.1, 1 a.a. ~ 527 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MSDTLTADVIGRRVEVNGEHATVRFAGVPPVAGPWLGVEWDNPERGKHDGSHEGTVYFK
CRHPTGGSFIRPNKVNFGTDFLTAIKNRYVLEDGPEEDRKEQIVTIGNKPVETIGFDSIM
KQQSQLSKLQEVSLRNCAVSCAGEKGGVAEACPNIRKVDLSKNLLSSWDEVIHIADQLRH
LEVLNVSENKLFPSGSLTGTLSVLKVLVLNQTGITWAEVLRVAGCPGLEELYLESNN
IFISERPTDVLQTVKLLDLSSNQLIDENQLYLIAHLPRLEQLILSDTGISLHFPDAGIG
CKTSMFPSLKYLVVNDNQISQWSFFNELEKLPSLRALSCLRNPLTKEDKEAETARLLIIA
SIGQLKTLNKCEILPEERRRAELDYRKAFGNEWKQAGGHKDPEKNRLSEEFHTAHPRYQF
LCLKYGAPEDWELKTQQPLMLKNQLLTLKIKYPHQLDQKVLEKQLPGSMTIQKVKGLLSR
LLKVPVSDLLLSYESPKKPGREIELENDLKSLLQFYVSVENGDCLLVRW


HRD, KCS, KCS1, pac2

Application WB

notes

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~59.30 kDa using TBCE transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of TBCE expression in transfected 293T cell line by TBCE rabbit polyclonal antibody.


Lane 1: TBCE transfected lysate(59.30 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TBCE transfected lysate.

Immunoprecipitation



Immunoprecipitation of TBCE transfected lysate using TBCE rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TBCE purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619