

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TSDP2-514

Anti TFDP2

Product type Primary antibodies

Immunogen TFDP2 (NP_006277.1, 1 a.a. ~ 386 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MIISTPQRLTSSGSVLIGSPYTPAPAMVTQTHIAEATGWVPGDRKRARKFIDSDFSESKR
SKKGDKNKGKGLRHFSMKVCEKVQRKGTTSYNEVADELVSEFTNSNNHLAADSQAYDQKNI
RRRVYDALNVLMAMNIISKEKKEIKWIGLPTNSAQECQNLEIEKQRRIERIKQKRAQLQE
LLLQQIAFKNLVQRNRQNEQQNQGPALNSTIQLPFIINTSRKTVIDCSISSDKFEYLF
NFDNTFEIHDDIEVLKRMGMSFGLESKGCSLEDLKLAKSLVPKALEGYITDISTGPSWLN
QGLLLNSTQSVSNLDLTTGATLPQSSVNQGLCLDAEVALATGQFLAPNSHQSSSAASHCS
ESRGETPCSFNDEDEEDDEEDSSSPE


Gene Alias DP2, Dp-2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.00 kDa using TFDP2 transfected 293T ystate.

Western Blotting (Transfected lysate)



Western Blot analysis of TFDP2 expression in transfected 293T cell line by TFDP2 rabbit polyclonal antibody.

Lane 1: TFDP2 transfected lysate(43.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619