

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TMEM177-305


Anti TMEM177

Product type	Primary antibodies
Immunogen	TMEM177 (AAH04404.1, 31 a.a. ~ 312 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	PVSYHLFPDPVVQWLYQYWPQGQPAPLPPQLQSLFQEVLQDIGVPSGHICYKPFTTFTFQP VSAGFPRLPAGAVVGIPASFLGDLVINTNHPVVIHGHTVDWRSPAGARLRASLTLSREAQ KFALAREVVCLESSTTAVHALLAPACLAGTWALGVGAKYTLGLHAGPMNLRAAFSLVAAV AGFVAYAFSQDSLTHAVESWLDLDRRTASLSAAYACGGVEFYEKLLSGNLALRSLLGKEGEK LYTPSGNIVPRHLFRIKHLPYTTRRDSVLQMWRGMLNPGRS*

Gene Alias MGC10993

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.02 kDa using TMEM177 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TMEM177 expression in transfected 293T cell line by TMEM177 rabbit polyclonal antibody.

Lane 1: TMEM177 transfected lysate(31.02 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619