

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-TP53-517**

# Anti TP53

**Product type** Primary antibodies

**Immunogen** TP53 (NP\_000537.2, 1 a.a. ~ 393 a.a) full-length human protein.

**Form** Liquid  
Protein A purified

**Volume** 100ul

**Specificity** Human

**Host** Rabbit

**Immunogen Sequence** MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQWFTEDPGP  
DEAPRMPEAAPRVAPAPAAPTPAAPAPAPSWPLSSSVPSQKTYQGSYGFRLGFLHSGTAK  
SVTCTYSPALNKMFCQLAKTCPVQLWVDSTPPPGRVVRAMAIYKQSQHMTEVVRRCPPHHE  
RCSDSDGLAPPQHILIRVEGNLRVEYLDDRNTFRHSVVVPYEPPEVGSDCTTIHYNMCNS  
SCMGGMNRRPILTIITLEDSSGNLLGRNSFEVRVCACPGRRRTEENLRKKGEPHHELP  
PGSTKRALPNNTSSSPQPKKKPLDGEYFTLQIRGRERFEMFRELNEALELKDAQAGKEPG  
GSRAHSSHLKSKKGQSTSRHKKLMFKTEGPDSD


**Gene Alias** FLJ92943, LFS1, TRP53, p53

**Application notes**

**WB**

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.23 kDa using TP53 transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of TP53 expression in transfected 293T cell line by TP53 rabbit polyclonal antibody.

Lane 1: TP53 transfected lysate(43.23 KDa).

Lane 2: Non-transfected lysate.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619