

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TP53I3-531

Anti TP53I3

Product type Primary antibodies

Immunogen TP53I3 (AAH00474.1, 1 a.a. ~ 332 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MLAVHFDKPGGPENLYVKEVAKPSPGEGEVLLKVAASALNRADLMQRQGQYDPPPGASNI
LGLASGHVAELGPCQGHWKIGDTAMALLPGGGQAQYVTVPEGLLMPIPEGLTLTQAAA
IPEAWLTAFQLLHLVGNVQAGDYVLIHAGLSGVGTAAIQLTRMAGAIPLVTAGSQKKLQM
AEKLGAAAGFNYKKEDFSEATLKFTKGAGVNLILDCIGGSYWEKNVNCALDGRWVLYGL
MGGGDINGPLFSKLLFKRGLITSLLRSRDNKYKQMLVNAFTEQILPHFSTEGPQRLLPV
LDRIYPVTEIQEAHKYMEANKNIGKIVLELPQ


Gene Alias PIG3

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~36.63 kDa using TP53I3 transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of TP53I3 expression in transfected 293T cell line by TP53I3 rabbit polyclonal antibody.


Lane 1: TP53I3 transfected lysate(36.63 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TP53I3 transfected lysate

Immunoprecipitation



Immunoprecipitation of TP53I3 transfected lysate using TP53I3 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TP53I3 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619