

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TRIM14-388

Anti TRIM14

Product type	Primary antibodies
Immunogen	TRIM14 (NP_150089.1, 1 a.a. ~ 442 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MAGAATGSRTTPGRSELVEGCGWRCPEHGDRVAELFCRRRCRCVCALCPVLGAHRGHPVGL
ALEAAVHVQKLSQECLKQLAIKKQQHIDNITQIEDATEKLKANAESSKTWLKKGKFTLRL
LLDEEEALAKKFIDKNTQLTLQVYREQADSCREQLDIMNDLSNRVWSISQEPDPVQRLQA
YTATEQEMQQQMSLGELCHPVPLSFEPVKSFFKGLVEAVESTLQTPDLIRLKESINCQLS
DPSSTKPGTLLKTSPSPERSLLLKYARTPTLDPDTMHARLRLSADRLTVRCGLLGSLGPV
PVLRFDALWQVLARDCFATGRHYWEVDVQEAGAGWWVGAAYASLRRRGASAAARLGCNRQ
SWCLKRYDLEYWAFHDGQRSRLRPRDDLRLGVFLDYEAGVLAFYDVTGGMSHLHTFRAT
FQEPLYPALRLWEGAISIPRLP

KIAA0129

Application notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.80 kDa using TRIM14 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TRIM14 expression in transfected 293T cell line by TRIM14 rabbit polyclonal antibody.

Lane 1: TRIM14 transfected lysate(49.80 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619