

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-TRIM22-157

Anti TRIM22

Product type	Primary antibodies
Immunogen	TRIM22 (NP_006065.2, 1 a.a. ~ 498 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MDFSVKVDIEKEVTCPICLELLTEPLSLDCGHSFCQACITAKIKESVIISRGESSCPVCQ TRFQPGNLRPNRHLANIVERVKEVKMSPQEGQKRDVCEHHGKKLQIFCKEDGKVICWVCE LSQEHQGHQTFRINEVVKECQEKLQVALQRLIKEDQEAEKLEDDIRQERTAWKNYIQIER QKILKGFNEMRVILDNEEQRELQKLEEGEVNVLNLAATDQLVQQRQDASTLISDLQRR LRGSSVEMLQDVIDVMKRSESWTLKKPKSVSKKLKSVFRVPDLSGMLQVLKELTDVQYYW VDVMLNPGSATSNAISVDQRQVKTVRTCTFKNSNPCDFSAFGVFGCQYFSSGKYYWEVD VSGKIAWILGVHSKISSLNKRKSSGFAFDPSVNYSKVYSRYRPQYGYWVIGLQNTCEYNA FEDSSSDPKVLTLMFMAVPPCRIGVFLDYEAGIVSFFNVTNHGALYKFSGCRFSRPAYP YFNPWNCLVPMTVCPSS
Gene Alias	GPSTAF50, RNF94, STAF50
Application notes	IP Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TRIM22 transfected lysate

Immunoprecipitation



Immunoprecipitation of TRIM22 transfected lysate using TRIM22 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TRIM22 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619