



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-TRIM72-164


Anti TRIM72

Product type	Primary antibodies
Immunogen	TRIM72 (AAH33211.1, 1 a.a. ~ 269 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSAAPGLLHQELSCPLCLQLFDAPVTAECGHSFCRACLGRVAGEPAADGTVLCPCCQAPT RPQALSTNLQLARLVEGLAQVPQGHCEEHLDPISYCEQDRALVCGVCASLGSHRGHRL PAAEAHARLKTQLPQQKLQLQEACMRKEKSVAVLEHQLVEVEETVRQFRGAVGEQLGKMR VFLAALEGLDREAERVGEAGVALRRELGLNSYLEQLRQMEKVLEEVADKPQTEFLMK YCLVTSRLQKILAESPPEARLDIQLPIIS

Gene Alias

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~29.50 kDa using TRIM72 transfected 293T lysate..
-------------------	--

Western Blotting (Transfected lysate)




Western Blot analysis of TRIM72 expression in transfected 293T cell line by TRIM72 rabbit polyclonal antibody.

Lane 1: TRIM72 transfected lysate(29.50 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TRIM72 transfected lysate.

Immunoprecipitation



Immunoprecipitation of TRIM72 transfected lysate using TRIM72 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TRIM72 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619