

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-UCK2-251

Anti UCK2

Product type Primary antibodies

Immunogen UCK2 (NP_036606.2, 1 a.a. ~ 261 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAGDSEQTLQNHQQPNGGEPFLIGVSGGTASGKSSVCAKIVQLLGQNEVDYRQKQVVILS
QDSFYRVLTSQKAKALKGQFNFDHPDAFDNELILKTLKEITEGKTVQIPVYDFVSHSRK
EETVTVPADVVLFEGLAFYSQEVRLDFQMKLFVDTDADTRLSRRVLRDISERGRDLEQ
ILSQYITFVKPAFEEFCLPTKKYADVIIIPRGADNLVAINLIVQHIQDILNGGPSKRQTNG
CLNGYTPSRKRQASESSSRPH

Gene Alias TSA903, UK, UMPK

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~29.30 kDa using UCK2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of UCK2 expression in transfected 293T cell line by UCK2 rabbit polyclonal antibody.

Lane 1: UCK2 transfected lysate(29.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619