



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-VTI1B-081


Anti VTI1B

Product type	Primary antibodies
Immunogen	VTI1B (AAH03142.1, 1 a.a. ~ 232 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MASSAASSEHFEKLHEIFRGLHEDLQGVPERLLGTAGTEEKKKLIRDFDEKQQEANETLA EMEEELRYAPPSFRNPMMSKLRNYRKDLAKLHREVRSTPLTATPGGRGGMKYGIYAVENE HMNRLQSQRAMLLQGTESLNRATQSIERSHRIATETDQIGSEIIEELGEQRDQLERTKSR LVNTSENLSKSRKILRSMRKYTTNKLLLSIIILLELAILGGLVYYKFFRSH

Gene Alias VTI1, VTI1-LIKE, VTI1L, VTI2

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~25.63 kDa using VTI1B transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of VTI1B expression in transfected 293T cell line by VTI1B rabbit polyclonal antibody.

Lane 1: VTI1B transfected lysate(25.63 KDa).

Lane 2: Non-transfected lysate.

Western Blotting (Cell lysate)



VT11B rabbit polyclonal antibody. Western Blot analysis of VT11B expression in Raw 264.7.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619