

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-027

Anti PRMT2

Product type	Primary antibodies
Immunogen	PRMT2 (NP_001526.2, 1 a.a. ~ 433 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MATSGDCPRSESQGEEPAECSEAGLLQEGVQPEEFVAIADYAATDETQLSFLRGEKILIL
 RQTTADWWWGERAGCCGYIPANHVGKHVDEYDPEDTWQDEEYFGSYGTLKLHLEMLADQP
 RTTKYHSVILQNKESLTDKVILDVGCCTGIISLFCAYARPRAVYAVEASEMAQHTGQLV
 LQNGFADIITVYQQKVEDVVLPEKVDVLVSEWMGTCLLFEFMIESILYARDAWLKEDGVI
 WPTMAALHLVPCSAADKYRSKVLFWDNAYEFNLSALKSLAVKEFFSKPKYNHILKPEDCL
 SEPCTILQLDMRTVQISDLETLRGELRFDIRKAGTLHGFTAWFSVHFQSLQEGQPPQVLS
 TGPFHPTTHWKQTLFMMDDPVPVHTGDVVTGSVVLQRNPVWRRHMSVALSWAVTSRQDPT
 SQKVGEKVFPIWR

Gene Alias HRMT1L1, MGC111373

Application notes **WB**
 Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.00 kDa using PRMT2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of PRMT2 expression in transfected 293T cell line by PRMT2 rabbit polyclonal antibody.

Lane 1: PRMT2 transfected lysate(49.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PRMT2 expression in Hela NE (Predicted Mw: 49.0420 kDa).

Western Blotting (Cell lysate)



PRMT2 rabbit polyclonal antibody. Western Blot analysis of PRMT2 expression in Hela NE.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619