

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-028

Anti DLAT

Product type	Primary antibodies
Immunogen	DLAT (NP_001922.2, 1 a.a. ~ 647 a.a) full-length human protein..
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MWRVCARRAQNVAPWAGLEARWTALQEVPGTPRVTSRSGPAPARRNSVTTGYGGVRALCG
WTPSSGATPRNRLLLQLLGSPGRRYYSLPPHQKVLPLSLSPTMQAGTIARWEKKEGDKIN
EGDLIAEVETDKATVGFESLEECYMAKILVAEGTRDVPIGAIICITVGKPEDIEAFKNYT
LDSSAATPQAAPAPTPAATASPPTPSAQAPGSSYPHMQVLLPALSPMTMTMGTVQRWEK
KVGEKLSEGDLLAEIETDKATIGFEVQEEGYLAKILVPEGTRDVPLGTPLCIIVEKEADI
SAFADYRPTEVTDLKPQVPPPTPPPVAAPVPTPQPLAPTPSAPCPATPAGPKGRVFSPL
AKKLAVEKGIDLTQVKGTGPDGRITKKDIDSFVPSKVAPAPAAVVPPTGPGMAPVPTGVF
TDIPISNIRRVIAQRLMQSKQTIPHYLLSIDVNMGEVLLVRKELNKILEGRSKISVNDFI
IKASALACKLVPEANSSWMDTVIRQNHVVDVSVAVSTPAGLITPIVFNAHIKGVETIAND
VVSLATKAREGKLQPHEFQGGTFTISNLGMFGIKNFSAIINPPQACILAIGASEDKLVPA
DNEKGFDVASMMSVTLSCDHRVVDGAVGAQWLAEFRKYLEKPITMLL


Gene Alias DLTA, PDC-E2, PDCE2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~69.00 kDa using DLAT transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of DLAT expression in transfected 293T cell line by DLAT rabbit polyclonal antibody.

Lane 1: DLAT transfected lysate(69.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DLAT expression in A-431 (Predicted Mw: 65.7810 kDa).

Western Blotting (Cell lysate)



DLAT rabbit polyclonal antibody. Western Blot analysis of DLAT expression in A-431.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619