

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-078

Anti CD1D


Product type	Primary antibodies
Immunogen	CD1D (AAH27926.1, 1 a.a. ~ 335 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MGCLLFLLWALLQAWGSAEVPQRLFPLRCLQISSFANSSWTRTDGLAWLGELQTHSWSN DSDTVRSLKPWSQGTFSQQWETLQHIFRVYRSSFTRDVKEFAKMLRSLYPLELQVSAGC EVHPGNASNNFFHVAFAQGKDILSFQGTSEPTQEAPLWVNLAIQVLNQDKWTRETVQWLL NGTCPQFVSGLLESGKSELKKQVKPKAWLSRGPSPGPGRLLLVCHVSGFYKPVWVKWMR GEQEQQGTQPGDILPNADETWYLRATLDVVAGEAAGLSCRVKHSSLEGQDIVLYWGGSYT SMGLIALAVLACLLFLLIVGFTSRFKRQTSYQGVL
---------------------------	--

Gene Alias CD1A, MGC34622, R3

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~36.96 kDa using CD1D transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CD1D expression in transfected 293T cell line by CD1D rabbit polyclonal antibody.

Lane 1: CD1D transfected lysate (36.96 kDa).

Lane 2: Non-transfected lysate.

1

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619