

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-098


Anti DDX50

[illegible]

Gene Alias	GU2, GUB, MGC3199, RH-II/GuB
-------------------	------------------------------

Application	WB
notes	Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~82.60 kDa using DDX50 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of DDX50 expression in transfected 293T cell line by DDX50 rabbit polyclonal antibody

Lane 1: DDX50 transfected lysate(82.60 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DDX50 expression in HeLa (Predicted Mw: 82.5650 kDa).


Western Blotting (Cell lysate)



DDX50 rabbit polyclonal antibody. Western Blot analysis of DDX50 expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of DDX50 expression in human kidney (Predicted Mw: 82.5650 kDa).

Western Blotting (Tissue lysate)



DDX50 rabbit polyclonal antibody. Western Blot analysis of DDX50 expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619