

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-109

Anti SPRY2

Product type Primary antibodies

Immunogen SPRY2 (NP_005833.1, 1 a.a. ~ 315 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MEARAQSGNGSQPLLQTPRDGGRQRGEPDPRDALTTQQVHVLSDQIRAIRNTNEYTEGPT
VVPRPGLKPAPRPSTQHKHERLHGLPEHRQPPRLQHSQVHSSARAPLSRSISTVSSGSRS
STRTSTSSSSSEQRLLGSSFSSGPVADGIIRVQPKSELKPGELKPLSKEDLGLHAYRCED
CGKCKCKECTYPRPLPSDWICDKQCLCSAQNVIDYGTCCVKGLFYHCSNDDDEDNCADN
PCSCSQSHCCTRWSAMGVMSLFLPCLWCYLPKAGCLKLCQGCYDRVNRPGCRCKNSNTVC
CKVPTVPPRNFEKPT

Gene Alias MGC23039, hSPRY2

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.70 kDa using SPRY2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SPRY2 expression in transfected 293T cell line by SPRY2 rabbit polyclonal antibody.

Lane 1: SPRY2 transfected lysate(34.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619