

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-133

Anti MPP5

Product type Primary antibodies

Immunogen MPP5 (NP_071919.2, 1 a.a. ~ 675 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MTTSHMNGHVTEESDSEVKNVDLASPEEHQKHREMAVDCPGDLGTRMMPIRRSAQLERIR
QQQEDMRRRREEEGKKQELDLNSSMRLKKLAQIPPKTGIDNPMFDTEEGIVLESPHYAVK
ILEIEDLFSSLKHIQHTLVDSQSQEDISLLLQLVQNKDFQNAFKIHNAITVHMNKASPPF
PLISNAQDLAQEVQTVLKPVHHKEGQELTALLNTPHIQALLAHDKVAEQEMQLEPITDE
RVYESIGQYGGETVKIVRIEKARDIPLGATVRNEMDSVVISRIVKGGAEEKSGLLEHGEDE
VLEINGIEIRGKDVNEVFDLLSDMHGTLTFVLIPSQQIKPPPAKETVIHVKAHFDYDPSD
DPYVPCRELGLSFQKGDILHVISQEDPNWWQAYREGDEDNQPLAGLVPGKSFQQQREAMK
QTIEEDKEPEKSGKLWCAKKNKKKRKKVLYNANKNDDYDNEEILTYEEMSLYHQPANRKR
PIILIGPQNCGQNELRQRLMNKEKDRFASAVPHTTSRRDQEVAGRQDYHFVSRQAFEDI
AAGKFIEHGEFEKNLYGTSIDSVRQVINSKICLLSLRTQSLKTLRNSDLKPYIIFIAPP
SQERLRALLAKEGKNPKPEELREIEKTREMEQNNGHYFDTAIVNSDLKAYQELLRLIN
KLDTEPQWVPSTWLR

Gene Alias FLJ12615, PALS1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~77.30 kDa using MPP5 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of MPP5 expression in transfected 293T cell line by MPP5 rabbit polyclonal antibody.

Lane 1: MPP5 transfected lysate(77.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619