

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-134

Anti SPRED1

Product type	Primary antibodies
Immunogen	SPRED1 (NP_689807.1, 1 a.a. ~ 444 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSEETATSDNDNSYARVRAVVMTRDDSSGGWLPLGGSGLSSTVFKVPHQEENGCAFFI RGERLRDKMVLVLECMKKDLIYNKVTPTFHHWKIDDKKFGLTFQSPADARAFDRGIRRAI EDISQGCPEKNEAEGADDLQANEEDSSSSLVKDHLFQQETVVTSEPYRSSNIRPSPFED LNARRVYMQSQANQITFGQPGLDIQSRSMYVQRQISKECGSLKSQNRVPLKSIRHVSFQ DEDEIVRINPRDILIRRYADYRHPDMWKNDLERDDADSSIQFSKPDSKKSDYLYSCGDET KLSSPKDSVVFKTQPSSLKIKKSKRRKEDGERSRCVYCQERFNHEENVRGKCQDAPDPIK RCIQVSCMLCAESMLYHCMSDSEGDFSDPCSCDTSDDKFCLRWLALVALSFIVPCMCCY VPLRMCHRCGEACGCCGGKHKAAG
Gene Alias	FLJ33903, NFLS
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~50.50 kDa using SPRED1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of SPRED1 expression in transfected 293T cell line by SPRED1 rabbit polyclonal antibody.

Lane 1: SPRED1 transfected lysate(50.50 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619