



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-158

Anti KATNA1


Product type	Primary antibodies
Immunogen	KATNA1 (AAH50428.1, 1 a.a. ~ 311 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MSLLMISENVKLAREYALLGNYDSAMVYYQGVLQDMNKYLYSVKDTYLQQKWQQVWQEIN VEAKHVKDIMKTLESFKLDSTPLKAAQHDLPASEGEVWSMPVPVERRRSPGPRKRQSSQY SDPKSHGNRPSTTVRVHRSSAQNVHNDRGKAVRCREKKEQNKGREEKGVLMMVGPPGTGKT LLAKAVATECKTTFFNVSSSTLTSKYRGESEKLVRLLFEMARFYSPATIFIDEIDSICSR RGTSEEHEASRRVKAELLVQMDGVGGTSENDDPKMMVMVLAATNFPWDIDEALRRRLEKR IYIPLPSGMRP
--------------------	--

Gene Alias

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.30 kDa using KATNA1 transfected 293T lysate.
-------------------	--

Western Blotting (Transfected lysate)



Western Blot analysis of KATNA1 expression in transfected 293T cell line by KATNA1 rabbit polyclonal antibody.

Lane 1: KATNA1 transfected lysate(35.30 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619