

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-190


Anti TNFAIP8L2

Product type	Primary antibodies
Immunogen	TNFAIP8L2 (AAH63014.1, 1 a.a. ~ 184 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MESFSSKSLALQAEKKLLSKMAGRSVAHLFIDETSSEVLDELRYVSKEYTHSRPQAQRVI KDLIKVAIKVAVLHRNGSFGPSELALATRFQRQKLRQGAMTALSFGGEVDFTFEAAVLAGLL TECRDVLLLELVEHHLTPKSHGRIRHVFDHFSDPGLLTALYGPDFTQHLGKICDGLRKLLD EGKL

Gene Alias FLJ23467, TIPE2

Application notes
WB
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~20.60 kDa using TNFAIP8L2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TNFAIP8L2 expression in transfected 293T cell line by TNFAIP8L2 rabbit polyclonal antibody.

Lane 1: TNFAIP8L2 transfected lysate(20.60 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619