



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-201

Anti SEC13

Product type	Primary antibodies
Immunogen	SEC13 (NP_899195.1, 1 a.a. ~ 322 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MVSVINTVDTSHEDMIHDAQMDYYGTRLATCSSDRSVKIFDVRNGGQILIADLRGHEGPV WQVAWAHPMYGNILASCSYDRKVIIWRENGTWEKSHEHAGHDSSVNSVCWAPHDYGLIL ACGSSDGAISLLTYTGEGQWEVKKINNAHTIGCNAVSWAPAVVPGSLIDHPGQKPNYIK RFASGGCDNLIKLVKEEEDGQWKEEQKLEAHSWVRDVAWAPSIGLPTSTIASCSQDGRV FIWTCDDASSNTWSPKLLHKFNDVVWHVSWSITANILAVSGGDNKVTLWKESVDGQWVCI SDVNKGQGSVSASVTEGQQNEQ
Gene Alias	d3S1231E, SEC13L1, SEC13R, npp-20
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.50 kDa using SEC13 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SEC13 expression in transfected 293T cell line by SEC13 rabbit polyclonal antibody.

Lane 1: SEC13 transfected lysate(35.50 kDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619