

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-204


Anti RTN3

Product type	Primary antibodies
Immunogen	RTN3 (NP_006045.1, 1 a.a. ~ 236 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAEPSAATQSHSISSSSFGAEPSPAGGGGSPGACPALGTKSCSSSCAVHDLIFWRDVKKT GFVFGTTLIMLLSLAAFSVISVVSYLILALLSVTISFRIYKSVIQAVQKSEEGHPFKAYL DVDITLSSEAFHNYMNAAMVHINRALKLIIRLFLVEDLVDSLKLAVFMWLMTYVGAVFNG ITLLILAEELLIFSVPIVYEKYKTQIDHYVGIARDQTKSIVEKIQAKLPGIAKKKA

Gene Alias ASYIP, HAP, NSPL2, NSPLI, RTN3-A1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~25.60 kDa using RTN3 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RTN3 expression in transfected 293T cell line by RTN3 rabbit polyclonal antibody.

Lane 1: RTN3 transfected lysate(25.60 KDa).

Lane 2: Non-transfected lysate.



COSMO BIO Co., LTD.

Inspiration for Life Science

Anti LC3 (Clone: 1703)

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619