

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-205


Anti ZCCHC11

Product type	Primary antibodies
Immunogen	ZCCHC11 (AAH48301.1, 1 a.a. ~ 309 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEESKTLKSENHEPKKNVICEESKAVQVIGNQTLKARNDKSVKEIENSSPNRNSSKKNKQ NDICIEKTEVKSCKVNAANLPGPKDLGLVLRDQSHCKAKKFPNSPVKAEKATISQAKSEK ATSLQAIAEKSPKSPNSVKA EKASSYQMKSEKVPSSPAEAEKGPSLLLKDMRQKTELQQI GKKIPSSFTSVDKVNIEAVGG EKCALQNSPRSQKQQTCTDNTGDSDDSASGIEDVSDDL KMKNDENKENSSEMDYLENATVIDESALTPEQRLGLKQAEERLERDHIFRLEKVYYVVL VIWGEMCVS

Gene Alias FKBP19, MGC54182

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.10 kDa using ZCCHC11 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of ZCCHC11 expression in transfected 293T cell line by ZCCHC11 rabbit polyclonal antibody.

Lane 1: ZCCHC11 transfected lysate(34.10 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619