

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-226

Anti MAN1C1

Product type	Primary antibodies
Immunogen	MAN1C1 (NP_065112.1, 1 a.a. ~ 630 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLMRKVPGFVPASPWGLRLPQKFLFLFLSGLVTLFCGALFLLPHSSRLKRLFLAPRTQQ PGLEVVAEIAGHAPAREQEPPPNPAPAAPAGEDDPSSWASPRRRKGGRLRRTRPTGPREE ATAARGNSIPASRPGDEGVPRFDNFNFRSRLRHPVLGTRADESQEPQSQVRAQREKIKE MMQFAWQSYKRYAMGKNELRPLTKDGYEGNMFGGLSGATVIDSLDTLYLMELKEEFQEAK AWVGESFHLNVSGEASLFEVNIRYIGLLSAFYLTGEEVFRIKAIRLGEKLLPAFNTPTG IPKGVVSFKSGNWWGATAGSSSILAEFGSLHLEFLHTELSGNQVFAEKVRNIRKVLRKI EKPFGLYPNFLSPVSGNWWQHHVSVGGGLGDSFYEYLIKSWLMSGKTDMEAKNMYEAELEA IETYLLNVSPGGLTYIAEWRGGILDHKMGHLACFSGGMIALGAEDAKEEKRAHYRELAQ ITKTCHESYARSDTKLGPEAFWFNSGREAVATQLSESYILRPEVVESYMYLWRQTHNPI YREWGWEVVLAEKYCRTEAGFSGIQDVYSSTPNHDNKQQSFFLAETLKLYLLFSEDDL LSLEDWVFNTTEAHPLPVNHSDSSGRAWGRH
Gene Alias	HMIC, MAN1A3, MAN1C, pp6318
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~70.90 kDa using MAN1C1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of MAN1C1 expression in transfected 293T cell line by MAN1C1 rabbit polyclonal antibody.

Lane 1: MAN1C1 transfected lysate(70.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619