



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-230

Anti ADSL

Product type	Primary antibodies
Immunogen	ADSL (NP_000017.1, 1 a.a. ~ 484 a.a) full-length human protein..
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAGGDHGSPPSYRSPLASRYASPEMCFVFSDRYKFRTWRQLWLWLAEAEQTLGLPITDE QIQEMKSNLENIDFKMAAEEEEKRLRHDVMAHVHTFGHCCPKAAGIIHLGATSCYVGDNTD LIILRNALDLLLPKLARVISRLADFAKERASLPTLGFTHFQPAQLTTVGKRCCLWIQDLC MDLQNLKRVRDDLRFGRGVKGTTGTQASFLQLFEGDDHKVEQLDKMVTEKAGFKRAFIITG QTYTRKVDIEVLSVLASLGASVHKICTDIRLLANLKEMEPEFEKQQIGSSAMPYKRNPMP SERCCSLARHMLTLVMDPLQTASVQWFERTLDDSANRRICLAEAFLTADTILNTLQNISE GLVVYPKVIERRIRQELPFMATENIIMAMVKAGGSRQDCHEKIRVLSQQAASVVKQEGGD NDLIERIQVDAYFSPHSQLDHLLDPSSFTGRASQQVQRFLIEEEVYPLLKPYESVMKVKA ELCL
Gene Alias	AMPS, ASASE, ASL
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~54.90 kDa using ADSL transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of ADSL expression in transfected 293T cell line by ADSL rabbit polyclonal antibody.

Lane 1: ADSL transfected lysate(54.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of ADSL expression in K-562 (Predicted Mw: 54.8890 kDa).

Western Blotting (Cell lysate)



ADSL rabbit polyclonal antibody. Western Blot analysis of ADSL expression in K-562.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619