

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-237

Anti PNLIP

Product type Primary antibodies

Immunogen PNLIP (AAH14309.1, 1 a.a. ~ 465 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MLPLWTLSLLLGA VAGKEVCYERLGCFSDDSPWSGITERPLHILPWSPKDVNTRFLLYTN
ENPNNFQEVAADSSSISGSNFKTNRKTRFIIHGFIIDKGEENWLANVCKNLFKVESVNCIC
VDWKGG SRTGYTQASQNI RIVGA EVAYFVEFLQSAFGYSPSNVHVIGHSLGAHAAGEAGR
RTNGTIGRITGLDPAEPCFQGTPELVRLDPSDAKFVDVIHTDGAPIVPNLGFGMSQVVGH
LDFFPNGGVEMPGCKKNILSQIVDIDGIWEGTRDFAACNHLRSYKYYTDSIVNPDGFAGF
PCASYNVFTANKCFPCPSGGCPQMGHYADRYPGKTNDVGQKFYLDTGDA SNFARWRYKVS
VTLSGKKVTGHILVSLFGNKGNSKQYEIFKGTLPDSTHSNEFSDVDVGD LQMVKFIWY
NNVINPTLPRVGASKIIVETNVGKQFNFCSPETVREEVLLTLTPC

Gene Alias PL

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~51.20 kDa using PNLIP transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of PNLIP expression in transfected 293T cell line by PNLIP rabbit polyclonal antibody.

Lane 1: PNLIP transfected lysate(51.20 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PNLIP expression in human pancreas (Predicted Mw: 51.1570 kDa).

Western Blotting (Tissue lysate)



PNLIP rabbit polyclonal antibody. Western Blot analysis of PNLIP expression in human pancreas.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619